# MORPHEUS M-E

ANAESTHESIA  
WORKSTATION  
FOR ADULTS,  
CHILDREN AND  
NEW BORN BABIES

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## MORPHEUS M

## MORPHEUS E

The Morpheus is an anaesthesia workstation suitable for treatment of adults, children and new-borns. It is provided in 2 versions: with electronic flowmeters group or with 5 rotameters keeping both the same characteristics of the system. The unit frame is made for 30% of light alloy and for 70% of plastic material. Thanks to the material used, the Morpheus has a total weight of less than 80 kgs.

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
The Morpheus LT is equipped with 3 wide and large drawers and a large work-shelf with LED lighting to facilitate operating rooms procedures.

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## ELECTRONIC VENTILATOR WITH 12" TFT COLOUR DISPLAY

The ventilator works by a microprocessor electronic control system with a simple and intuitive interface, which allows the anesthetist to operate in a clear and safe way. Unique control by encoder for the regulation and management of all respiratory parameters. This system incorporates the automatic compensation of dead space, the leakages calculation and the automatic calibration of  $\text{FiO}_2$  transducer.

### Operative modes - Modi operativi:

MANUAL,  
VC-VAC,  
APCV,  
PCV-TV,  
SPONT,  
SIMV+PS (volume),  
PSV,  
APNOEA BACK-UP.

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### Ventilation curves:

Pressure  
Flow  
Volume  
Loop (volume/pressure, volume/flow)

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Paw  
PEEP  
Vte  
ExpMV  
Frequenza  
I:E  
 $\text{FIO}_2$   
Pmax  
Pmedia  
Pausa (%)  
Vti  
InspFlow  
ExpFlow  
Tinsp  
Tesp  
Tpausa  
Cs  
Ri

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Pressione  
Flusso  
Volume  
Loop (volume/pressione, volume/flusso)

### Electronic Flowmeters Group

Thanks to its ergonomic and compactness the electronic flowmeters group is user friendly and intuitive.

It allows the electronic selection of medical gas (O<sub>2</sub>- N<sub>2</sub>O-AIR) and it is also predisposed for use of Xenon gas.

It is equipped with safety scaled flowmeter, gauges for medical gas from central system, visual and acoustic alarms for gas supply hose disconnection and gas supply failure, back-up battery for 120 minutes operation.

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# anaesthesia workstations

### Mechanical Flowmeters Group

Thanks to its ergonomics and compactness the 5 rotameters flowmeters group is user friendly and intuitive.

It allows the electronic selection of medical gas (N<sub>2</sub>O and AIR), it is equipped with gauges for medical gas from central system, and visual and acoustic alarms for gas supply failure and gas supply hose disconnection.

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Medical gas connection from main gas supply system and cylinders, auxiliary outlets.

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## PATIENT CIRCUIT GROUP


The Morpheus breathing system is an absolute novelty. Its design philosophy is the result of almost 40 years experience in anaesthesia field. The system is compact, lightweight for manual removal and fully autoclavable. The patient group is supplied with an automatic heating system avoiding humidity with possibility to use IS PAC disposable CO<sub>2</sub> absorbers of 0.5 kgs or reusable and rechargeable CO<sub>2</sub> absorbers (TOP SPECIAL of 1.0 kgs.). Both absorbers are with quick-release enabling replacement also during anaesthesia. The group is supplied with:

- a bellows that can be manually inspected
- fresh gas collection bag
- FiO<sub>2</sub> cell
- MICRO FLOW<sup>®</sup> sensor
- Safety valve for collecting ambient air in case of lack of fresh gas

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## SIARE FLOW

The respiratory dynamics analysis of Morpheus is made through the exclusive MICRO FLOW<sup>®</sup> sensor with operating principle by magnetic disturbance. This reusable and sterilisable flow sensor is placed directly inside the expiratory valve and thanks to its advanced technology it can analyse volumes below 10 ml.

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The Morpheus is equipped with an advanced manual ventilation system that can be switched on and off automatically by the ventilator; the operator can select the APL system or the AUX system with traditional To-and-Fro valve and O<sub>2</sub> by pass safety system.

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### MORPHEUS M

#### OM3 S5 MORPHEUS M

Anaesthesia Workstation including:

- Mechanical gas mixing system at 3 gas and 5 rotameters (O<sub>2</sub>, N<sub>2</sub>O and Air).
- Electronic ventilator.
- Patient group for ventilation in open, semi-closed, closed at low flow system.
- Support for the rapid connection of two vaporizers.
- TOP SPECIAL CO<sub>2</sub> absorber.

### MORPHEUS E

#### OM3 SE MORPHEUS E

Anaesthesia Workstation including:

- Electronic gas mixing system at 3 gas (O<sub>2</sub>, N<sub>2</sub>O, Air).
- Electronic ventilator.
- Patient group for ventilation in open, semi-closed, closed at low flow system.
- Support for the rapid connection of two vaporizers.
- TOP SPECIAL CO<sub>2</sub> absorber.

**Code 100150/IDF**  
Isoflurane  
Vaporizer

**Code 100148/IDF**  
Fluo-Halotane  
Vaporizer

**Code 100145/IDF**  
Sevoflurane  
Vaporizer

**Code 100151/IDF**  
Desflurane  
Vaporizer

**Code 06751**  
Stainless steel  
patient circuit  
supporting arm

**Code 30005/10**  
IS PAC  
0,5 Kg CO<sub>2</sub> absorber disposable  
BOX of 10 pcs

**Code 30900**  
TOP SPECIAL  
1 Kg CO<sub>2</sub> absorber rechargeable

**Code A56.NEPT.15**  
NEPTUNE 15  
Multi-Parameter  
Patient Monitor  
with 15" TFT display

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REG. N. 327  
UNI EN ISO 9001-2000  
UNI CEI EN ISO 13485-2002

CE 0476

SIARE is not only research, designing and manufacturing but with an efficient and expertise technical service around the world grants the maximum efficiency of its products.

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