


# Siare



## NEPTUNE IPER

VITAL SIGNS MONITOR FOR HYPERBARIC OXYGEN THERAPY (HBOT)

- Matchless performance up to a depth of 60 mH<sub>2</sub>O
- Monitored parameters: ECG 3/5 lead – RESP – NIBP – TEMP (2 channels) SPO<sub>2</sub> – ARR ANALYSIS – ST ANALYSIS
- NIBP technology based on a pneumatic valve supplied with compressed AIR.
- Suitable for Adult, Paediatric and Neonatal monitoring
- 12" TFT LED COLOUR MONITOR with HIGH RESOLUTION.
- Small dimensions and light weight (3,5 Kg)
- Special hermetic battery with autonomy from 2 to 10 hours



**siare**  
ENGINEERING  
INTERNATIONAL GROUP s.r.l.

**20 years**  
**of proven experience**  
**in Hyperbaric Medicine**